

# Using metronomes & microphones to model a pulsar timing array





### Metronome Timing Array

microphone motion causes pulses to arrive ahead or behind schedule, correlated across metronomes (pulse perturber)

metronome (pulse emitter)

(pulse detector)

microphone

metronome (pulse emitter)

Key quantities to calculate


pulse period T<sub>P</sub>

pulse shape (profile)

measured times of arrival (TOAs) expected TOAs

residuals = measured - expected

correlations between pairs of residuals



#### Some numbers to keep in mind

```
v_{\text{sound}} = 340 \text{ m/s (in air)}
```

amplitude  $\approx 10$  cm

amplitude /  $v_{sound} = 3 \times 10^{-4} \text{ sec}$ 


120 bpm:  $T_p = 0.5 \text{ sec}$ 

200 bpm:  $T_p = 0.3 \text{ sec}$ 

### PTAdemo0GUI.py


m200a.txt m120b.txt

# Data folding to determine pulse period and pulse profile





- break data into chunks of duration T
- fold (stack) chunks, then add together
- get maximum sum when T is the period of the signal (noise cancels, signal combines coherently)

#### Metronome pulse profiles



# Estimating TOAs by correlating data with the pulse profile






### PTAdemo2GUI.py

m200a I 20b0.txt

### Expected residuals for circular microphone motion



### PTAdemo2GUI.py

m200a | 20b | 80.txt m200a | 20b 90.txt m200a | 20b 45.txt

m200a I 20b I 35.txt

### Metronome timing array illustrates several techniques from pulsar timing

- I. Estimate pulse period and pulse shape (profile) by folding data
- 2. Estimate pulse times of arrival (TOAs) by correlating data with pulse profile
- 3. Calculate timing residuals by subtracting expected TOAs (based on a model) from the measured TOAs
- 4. Improve estimate of pulse period by removing a linear trend from residuals
- 5. Timing residuals for a pair of pulsars are correlated as a function of their angular separation